




Stockholm
University


Enantioselective cyclopropane assembly

Key Organics is very pleased to exclusively offer a new reagent developed by Dr Abraham Mendoza and colleagues at Stockholm University as featured in the highly rated publication, *General Cyclopropane Assembly via Enantioselective Redox-Active Carbene Transfer to Aliphatic Olefins*

Marc Montesinos-Magraner, Matteo Costantini, Rodrigo Ramirez-Contreras, Michael E. Muratore, Magnus J. Johansson, Abraham Mendoza, *Angewandte Chemie*, Volume 131, Issue 18, Pages 5991-5996.



Within the paper, Dr Mendoza and colleagues introduce a practical and versatile diazocarbene, and demonstrate its performance in the first unified asymmetric synthesis of functionalized cyclopropanes. NHPI-DA(13a) is available from Key Organics as SO-3001.



EXCEPTIONAL SERVICE IN THE USA – EXCEPTIONALLY CLOSE BY

Our Bedford, MA facility is only 45 minutes' drive from the pharma, biotech, and academic hubs of Cambridge, Boston and Worcester. This means that we can get any of our 10,000 in-stock compounds over to your project very quickly. We offer same time zone responsiveness, so we can assist with synthesis and inventory questions on the very same day. The west coast hive of activity is covered too! Same day quotes, information and quick shipping including overnight service are available at reasonable rates.



Give us a call or email us:

Toll Free

855 808-2700

or 781 280-5000

E: enquiries@keyorganics.net

Key Organics Inc., Suite 100, 55 North Road, Bedford, MA 01730


Our growing expertise in Stable labelled API's and Intermediates

One currently popular area for our clients is the synthesis of cold (non-radioactive) labelled reference standards. These materials are of use as tracers in the modeling of a number of chemical and biochemical systems, for example, in tracking the passage of an isotope, or an atom with a variation, through a reaction, metabolic pathway, or cell. The reactant is 'labeled' by replacing specific atoms by their isotope. The position of the isotopes in the products is measured to determine the sequence the isotopic atom followed in the reaction or the cell's metabolic pathway.

[²H₄]-Linagliptin Case Study

We have developed considerable experience in the synthesis of ²H- (deuterium) and ¹³C-labelled materials. Whilst it is often possible to identify a suitably labelled commercially available building block with which to replicate the known non-labelled chemistry, there are instances where we need to devise an alternative synthesis route. Two such examples are to facilitate the late stage incorporation of the label, minimising the quantity – and cost – of material required, or where an appropriate building block is not available from the commercial pool.


In one particular example, a customer required gram scale quantities of [²H₄]-linagliptin. Linagliptin (CAS 668270-12-0) is a DPP-4 inhibitor developed by Boehringer Ingelheim for treatment of type II diabetes.¹ The availability of the late stage bromoimidazole building block **1** suggested placing the labelled atoms in the aminopiperidine component **2** (Scheme 3).




Scheme 3: [²H₄]-Linagliptin retrosynthesis

Adaption of a known method for the synthesis of enantiomerically pure 3-amino cyclic amines afforded the required building block (Scheme 4).² Commercial D-glutamic acid was esterified and N-Boc protected. Reduction of the two ester functions to the corresponding alcohols proceeded cleanly with sodium borohydride, and no loss of the Boc group was observed. Utilising the corresponding [²H₄]-reductant afforded the chiral [²H₄]-diol, which was converted to the bis-mesylate and cyclised to (R)-3-(N-Boc-amino)-2,2,4,4-[²H₄]piperidine on treatment with ammonia.

Scheme 4: (R)-3-(N-Boc-amino)-2,2,4,4-[²H₄]piperidine synthesis



With gram quantities of the labelled building block in hand, the synthesis of [²H₄]-linagliptin was completed via reaction of the customer supplied bromoimidazole building block **1** with the aminopiperidine component **2**, followed by Boc-deprotection (Scheme 5).¹



References:

- M. Eckhardt, E. Langkopt, M. Mark, M. Tadayyon, L. Thomas, H. Nar, W. Pfrengle, B. Guth, R. Lotz, P. Sieger, H. Fuchs and F. Himmelsbach, *J. Med. Chem.*, 2007, 50, 6450-6453
- S.-H. Moon and S. Lee, *Syn. Commun.*, 1998, 28, 3919-3926

Key Organics Ltd. & NTZ Lab Ltd. entered into advanced *in vivo* studies in relevant Parkinson's and Alzheimer's disease models

A new collection of structurally optimized compounds addressing the central nervous system (CNS) diseases, which have been recently entered into a development and marketing agreement between Key Organics Ltd. and NTZ Lab Ltd., is now available in the BIONET collection. The aim of the collaboration is to explore global markets for new generation of structurally optimized compounds addressing CNS diseases. This screening collection with experimentally-determined bioactivity, ADMET, bioavailability, and toxicity will address the need of new, robust, and multi-potent small molecules for the treatment and diagnosis of CNS diseases, such as Parkinson's disease, Alzheimer's disease, dementia and/or other neurodegenerative diseases.

The initial set of nine target indazole-5-carboxamides (subclasses I–III) is part of large number CNS active compounds. Importantly, the newly discovered molecules are virtually and experimentally evaluated as multi-target acting monoamine oxidase (MAO)/histone deacetylase 2 (HDAC2) inhibitors by using a well validated, combined X-ray/modeling technology platform (Fig.1). In addition, the neuroprotective effects on TH-positive dopaminergic neurons and the induction of the neurite network outgrowth of the most promising compound NTZ-1441 have been investigated in a Parkinson's disease model (Fig. 2). These effects are associated with a good BBB penetration of all presented compounds that was confirmed in several *in vitro* assays.*

Due to the excellent ADMET profile combined with experimental confirmed "Proof-of-Principle" biological activity, selected representatives of this next generation multi-target compounds are now advanced for further *in vivo* studies in relevant Parkinson's and Alzheimer's disease models.

Fig. 1: Multipotency Screening (A) and estimated binding affinity (K_i HYDE ranges) at hMAO-B (B) and hHDAC2 (C).





Fig. 2: Effect of NTZ-1441 on dopaminergic neurons and neurite network after 48 hrs injury with MPP+.



* $p < 0.05$ vs. MPP+ condition. Data are expressed as the mean \pm SEM ($n = 4-6$) and were analyzed using one-way ANOVA followed by Dunnett's test.

References:

*Carboxamides vs. methanimines: Crystal structures, binding interactions, photophysical studies, and biological evaluation of (indazole-5-yl)methanimines as monoamine oxidase B and acetylcholinesterase inhibitors. *European Journal of Medicinal Chemistry*, Volume 179, 1 October 2019, Pages 404-422. The article is available online at: <https://authors.elsevier.com/a/1ZJ3xoqIQORRP>

We will be exhibiting and/or attending the following exhibitions and conferences in the coming months:

| | | |
|----------------|---|----------------|
| August 21st | 27th Annual GP2A Medicinal Chemistry Conference | Nottingham, UK |
| September 8th | 20th RSC / SCI Cambridge MedChem symposium | Cambridge, UK |
| September 16th | Discovery On Target | Boston, US |

For more information, please visit: <https://www.keyorganics.net/about/exhibition-conference-attendance/>

Staff Interview Lorna Bankole, Analytical Chemist



Q: Please tell us a bit about yourself?

A: I have a BSc Honours degree in Applied Biochemistry from the University of Paisley with an MSc degree in Pharmaceutical Analysis from the University of Strathclyde along with 22 years of experience (& still counting) in the pharmaceutical/chemical industry. I have a good background in chromatography (HPLC/GC) & spectroscopy (UV, IR, NMR, mass spec). My hobbies/interests include

roller skating, watching TV (quiz shows & hospital dramas), listening to music & playing computer games. I am particularly passionate about Morris Minor cars & am pleased to say that after around 30 years of interest in the Morris Minor I am now the proud owner of a 1959 Morris Minor 4 door saloon (as seen in the photo above).

Q: What is your role within Key Organics?

A: I am an analytical chemist & my main role is to test samples (either internal or external) by LCMS, NMR, GC.

Q: What do you enjoy about working at Key Organics?

A: The people I work with are great & I get a lot of interaction & involvement within my role in addition to the varied aspect of the role.

Q: What do you think is Key Organics' greatest strength?

A: We strive to maintain on time delivery, providing a good quality service for both internal & external customers.

Key Organics
Chemistry | Innovation | Quality

For more information, please contact us at:

Key Organics Ltd.,
Highfield Road Industrial Estate,
Camelford,
Cornwall PL32 9RA,
UK

T: +44 (0)1840 212137

F: +44 (0)1840 213712

E: enquiries@keyorganics.net



Key Organics Inc.,
Suite 100,
55 North Road,
Bedford,
MA 01730

Toll Free No: 855 808-2700

Or: 781 280-5000

E: enquiries@keyorganics.net

www.keyorganics.net